

LIFEPAK® 12 defibrillator/monitor Series OPERATOR'S CHECKLIST

Unit Serial No.: _____
Location: _____

This is a suggested checklist to inspect and test this device. Inspection and testing of this device on a daily basis is recommended. You may also consult JAMA, August 22/29, 1990, Vol. 264, No. 8, Table 2 for the Defibrillator Working Group's manual defibrillator checklist.

This form may be reproduced.



Instruction	Recommended Corrective Action	Date							
		Initials							

Insert a ✓ in the box after completing each instruction.

1 Inspect physical condition for:

Foreign substances

Clean the device.

Damage or cracks

Contact qualified service technician.

2 Inspect Power Source for:

Broken, loose, or worn battery pins

Contact qualified service technician.

Damaged or leaking battery

Discard/recycle battery

Two fully charged batteries installed

Replace battery.

Fully charged spare batteries available

Obtain battery.

3 Check therapy and ECG electrodes for:

Expiration date

Replace if expired.

Spare electrodes available

Obtain spare electrodes.

4 Examine accessory cables for cracking, damage, broken or bent parts or pins, and paddle surfaces for pitting.

Replace damaged or broken parts.

5 With the battery installed and the device disconnected from ac power, press ON and look for:

Self-test messages

If absent, contact qualified service technician.

Momentary illumination of each LED and all LCD segments

If absent, contact qualified service technician.

LOW BATTERY or REPLACE BATTERY messages

Replace the battery immediately.

Service indicator message

Contact qualified service technician.

Perform User Test (QUIK-COMBO therapy cable only)

If test failed, repeat. If failed twice, contact qualified service technician.

6 Check ECG printer for:

Adequate paper supply.

Replace if necessary.

Ability to print.

If not working, contact qualified service technician.
