

Sticking Points:



COMMON CONCERNS SEEN BY EMSQA

Paramedics:

This is your opportunity
to lead up the chain



Initiate the conversation!

Teach your officers (CPR Bosses) how to use Handtevy

Make sure you bring your Sonim Phone in to cardiac arrests and give it to them

Teach them how to use the CPR timer, enter medication dosages, procedures, events. Teach them how to find equipment sizes for you.

Equipment	Instruction
Bag Resuscitator with PEEP Valve	Child (with manometer)
Laryngoscope Blade	McGrath 2 - Mac 2 - Mi 2
ETT (Have smaller and larger tubes ready)	4.5 - 5.0 Cuffed
Stylet	10 French
Suction Catheter	10 French
ETT @ Lip	15 - 16 cm
OPA (Teeth to Angle Jaw)	70 mm (Size 2)
NPA (Nostril to Earlobe)	24 French
LMA	Size 2 or 2.5
F-gel Supraglottic Airway	Size 2
King Laryngeal Tube	Size 2.5