

# LIFEPAK<sup>®</sup> 15 MONITOR/DEFIBRILLATOR

Works like you work.™

## Operationally Innovative, Technologically Advanced

The right tool for the job is critical, particularly for rescue personnel. While good technology can provide needed information, great technology works so intuitively that you never give it a second thought. LIFEPAK products are built with operational innovations that free you to focus entirely on your lifesaving mission—not your equipment.

Medical emergencies can happen anywhere—indoors or outside, in any type of weather, at any time of day—the LIFEPAK 15 monitor/defibrillator is intelligently designed with new innovations to support a multitude of emergency situations. No matter the location or environmental conditions, our commitment to operational innovation has yielded the best possible equipment for saving lives.

### What makes something operationally innovative?

LIFEPAK products are a lifeline for your patients—and a vital source of information for you. We understand these basic fundamentals, and design products to improve your speed, ease and effectiveness in saving lives on the job. We focused on four key essentials we know are important to you.

#### Viewability

Because patient treatment depends on the information displayed on the monitor screen, the 15 is designed to provide maximum viewability. You need the ability to read information in all kinds of situations and environments, and to see as much patient data as possible.

The LIFEPAK 15 monitor/defibrillator is equipped to meet your needs with:

- A large 8.4 inch diagonal display (6.7 W X 5 H)
- SunVue™ feature makes it easy to see in bright sunlight—the touch of a button eliminates washout
- Logical user interface with parameters arranged down one side and color matched to waveforms
- Largest combination of font and screen text of any monitor/defibrillator on the market



“When you’re fighting a fire you need a fire hose and a fire engine. When you go to a cardiac arrest you need the LIFEPAK.”

Richard Bylicki, Commissioner,  
South Farmingdale Fire Department,  
New York

## Power you can count on

Lives count on you as well as your equipment. The 15 is powered by the latest technology to help you know it's ready to go throughout your full shift.

- Two long-lasting Lithium-ion batteries have a total operational time of up to six hours.
- No need to wait for full recharge before using the device. Charge when you can, as you can.
- Smart battery technology displays the remaining battery life to help you manage inventory.

## A platform to last

We realize the LIFEPAK 15 monitor/defibrillator is a large expenditure for your service. We also know you want to keep current with the latest technologies. That's why we've made the 15 upgradeable so as standards evolve you can add new capabilities, over time, and extend the value of your investment.

- Increased processing speed and memory enable you to add new features and monitoring parameters; the platform is ready for anticipated Guidelines 2010 changes.
- Compatible with all leading electronic patient records systems to ensure vital patient information can be transferred from the 15 to your records system.
- Designed so information easily flows where it's needed for evaluation, summarizing and reporting.

## Industry Service Leader

The LIFEPAK 15 monitor/defibrillator is backed by the largest and most experienced service team in the industry. Our service representatives have an average tenure of over 12 years and log an average of 2,400 field hours each year. With comprehensive product coverage and 24-hour response time in most areas, you can count on the Physio-Control service team to support you when and where you need it.

## Design details that make a difference

The 15 is a hands-on tool so it must be durable, comfortable to work with, and easy to maintain.

- Buttons are glossy black on a matte black background and are easy to find, use and clean.
- Front-based connectors enable you to monitor with confidence, knowing you can see if a connection is properly connected.
- A larger SPEED DIAL makes it easy to make your selection.
- A wide ergonomic handle makes the 15 comfortable to carry and pass off to your EMS and hospital colleagues as needed.
- A new "power-off" button on the 15 is designed to reduce the chance of accidental power loss when moving the unit with a patient.

### Things for LIFEPAK 12 defibrillator/monitor users to consider:

Budget cycles and resources can sometimes constrain your goal of replacing your entire monitor fleet at once. The 15 was purposely designed to help reduce training and switching costs.

- A similar footprint to the LIFEPAK 12 defibrillator/monitor means there is already a spot earmarked for the unit in your rig—no need to reconfigure.
- Similar functionality and familiar user interface allows for easy transition to the 15 for your medics, minimizing training costs.

All statements and comparative information in this brochure are valid as of February 2009.

For further information, please contact your sales representative at 800.442.1142 or visit our website at [www.physio-control.com](http://www.physio-control.com)

**PHYSIO  
CONTROL**

#### HEADQUARTERS / MANUFACTURING

Physio-Control, Inc.  
11811 Willows Road NE  
P. O. Box 97006  
Redmond, WA 98073-9706 USA  
Tel 425 867 4000  
Fax 425 867 4121  
[www.physio-control.com](http://www.physio-control.com)

#### SALES OFFICES

**Europe**  
Medtronic International  
Trading Sàrl  
Case postale 84  
Route du Molliau 31  
CH-1131 Tolochenaz  
[www.physio-control.com](http://www.physio-control.com)  
Tél +41 (0)21 802 70 00  
Fax +41 (0)21 802 79 00

**Physio-Control UK**  
Medtronic Limited  
Building 5  
Croxley Green Business Park  
Watford, Hertfordshire  
WD18 8WW  
Tel 44 1923 212 213  
Fax 44 1923 241 004